

Product information

Product name:	anti- rubellin pab ELISA-Stick with test reagents
Product number:	ER-4
Package size:	5 tests
General Description:	Rubellin is a toxic anthraquinone derivative produced by <i>Ramularia collo-cygni</i> , the fungus causing <i>Ramularia</i> leaf spot disease on barley and some other cereals. <i>Ramularia collo-cygni</i> is now recognized as an important pathogen of barley in Northern Europe and New Zealand. The fungus produces a number of anthraquinone toxins (rubellins), which act as host nonspecific toxins with photodynamic activity. These toxins induce lipid peroxidation and are possibly the cause of the chlorosis and necrosis observed in leaves infected with <i>R. collo-cygni</i> . The fact that the fungus can remain latent in barley plants until flowering, coupled with its very slow growth in vitro, makes it difficult to detect in crops. As a result, the epidemiology of this pathogen remains poorly understood.
Principle:	The <i>Rubellin detection reagent set</i> is sensitive and easy to handle test for the detection of <i>Rubellin</i> . The test can be carried out without sample preparation or pre-incubation. Since the test does not require any labour equipment, it is especially suitable for analysis on the spot.
Sensitivity:	Detection limit 1 ng
Sample size:	max. 1 ml
Assay sample type:	any sample
Detection time:	1.5 h

All products sold by AntiProt are intended for research use only unless otherwise indicated. This product is not intended for diagnostic or drug purposes, or for use in humans.